


LC/MSMS Methods for Quantitative Determination of Mycophenolic Acid in Human Plasma

Mycophenolic acid (F1) is an immunosuppressive agent used in therapy following solid organ transplantation. The active drug component is marketed as the sodium salt (Myfortic, Novartis) or as an ester (Cellcept, Roche) that is metabolized *in vivo* to mycophenolic acid.

F1. Chemical Structure of Mycophenolic Acid



Liquid chromatography-based methods using either UV or mass spectroscopic detection is the standard for determination of mycophenolitic acid (MPA) in biological samples. Patel et. al.[1] reported using 100 mL sample extracted by SPE method analyzed by reverse phase chromatography with either UV detection (18 minute cycle time) or MS/MS detection (7 minute cycle time) with a quantification range of 1-1000 mg/mL. Analysis of plasma and plasma ultrafiltrate used LC-MS/MS with detection range of 0.5-1000 mg/mL range[2]. LC-UV method used 50 mL sample size with protein precipitation, and chromatographic separation on a monolithic column with isocratic elution gave linear response for MPA at 1-40 mg/mL range[3]. Fluorescence detection was used to analyze 100 mL plasma sample extracted by protein precipitation with methanol or direct injection of plasma ultrafiltrate with quantification range of 5-1000 mg/mL[4]. A method that included analysis of urine used reverse phase chromatography with ion-pairing by UV detection with quantification range of 5-400 mg/mL[5].

Stability of mycophenolate acyl glucuronides was reported in biological samples prepared for analysis by protein precipitation procedure followed by reverse phase HPLC separation with UV detection. (6)



T1. Inter-run (n=14) Accuracy and Precision for Standard Calibrators

Concentration (ng/mL)	SC 50 50.0 ng/mL	SC 100 100 ng/mL	SC 250 250 ng/mL	SC 1000 1000 ng/mL	SC 2500 2500 ng/mL	SC 5000 5000 ng/mL	SC 7500 7500 ng/mL	SC10000 10000 ng/mL
Mean	50.2	99.2	248	1000	2570	5120	7220	9900
%CV	3.8	4.8	2.7	3.6	3.5	3.2	3	2.9
%Bias	0.4	-0.8	-0.8	0	2.8	2.4	-3.7	-1
n	14	14	14	14	14	14	14	13



T2. Inter-run (n=18) Accuracy and Precision for Quality Control Samples

	Curve Number	QC 50 50.0 ng/mL	QC 150 150 ng/mL	QC 500 500 ng/mL	QC 8000 8000 ng/mL
Intrarun Mean	1	44.2	142	508	7360
Intrarun %CV		6.3	4.7	3.0	0.9
Intrarun %Accuracy		88.4	94.7	101.6	92.0
n		6	6	6	6
Intrarun Mean	3	54.4	153	506	7660
Intrarun %CV		4.1	1.5	1.5	2.7
Intrarun %Accuracy		108.8	102.0	101.2	95.7
n		6	6	6	6
Intrarun Mean	4	53.3	154	487	7750
Intrarun %CV		5.8	2.6	2.0	8.2
Intrarun %Accuracy		106.6	102.7	97.4	96.9
n		6	6	6	6
Inter-run Mean		50.6	150	500	7590
Inter-run %CV		10.6	4.8	2.9	5.3
Inter-run %Accuracy		101.2	100.0	100.0	94.9
n		18	18	18	18

F2. Extracted Blank



F3. LLOQ Sample (50 ng/mL)





POWERFUL RESULTS
You invested in discovery and innovation, now power your drug development with our 3 laboratories.

**45 LC/MSMS
GLP Bioanalysis
Support for Preclinical and
Clinical Pharmacokinetics**

100 analytical scientists, producing results of the highest quality and regulatory reputation with the latest in automation, data acquisition, interpretation and reporting

BASi
West Lafayette, IN • McMinnville, OR • Warwick, UK
www.bioanalytical.com • basi@bioanalytical.com • 800.845.4246

F4. ULOQ Sample (10000 ng/mL)



BASi has developed and validated an LC/MSMS method for analysis of mycophenolic acid (immunosuppressant drug) in human plasma. The method uses protein precipitation for quick sample preparation. Processing a sample volume of only 0.1 mL plasma yields a detection limit of 50 ng/mL. The method for analysis of mycophenolate plasma was shown to be robust with standard calibrator accuracy (CV) well within acceptance guidelines of $\pm 15\%$ over the course of seven independent sequences (T1). The between-run accuracy and precision of quality control samples prepared at four different levels [LLOQ (QC-50), low-QC (QC-150), mid-QC (QC-500) and high-QC (QC-8000)] run in three sequences with $n=6$ replicates in each sequence also showed excellent method performance (T2). Representative chromatograms are shown in F2 thru F4. Ruggedness tests of freeze/thaw stability, heat treatment stability (for deactivation of HIV virus), short-term room temperature stability and autosampler stability were performed with acceptable accuracy and precision. No matrix interferences were detected testing six different lots of plasma.

This analytical method for determination of MPA is a reliable and convenient procedure that meets the criteria for application in routine clinical drug monitoring and pharmacokinetic studies.

Reference

1. C. G. Patel, A. E. Mendonza, F. Akhlaghi, O. Majib, A. K. Trull, T. Lee, D. W. Holt *J Chrom B*, 813 (2004) 287-294.
2. F. Streit, M. Shipkova, V. W. Armstrong, M. Oellerich, *Clin Chem* 50:1 (2004) 152-159.
3. D. Indjova, L. Kassabova, D. Svinarov, *J Chrom B*, 817 (2005) 327-330.
4. J. Shen, Z. Jiao, Y. Yu, M. Zhang, M. Zhong, *J Chrom B*, 817 (2005) 207-213.
5. W. Yau, A. Vathsala, H. Lou, E. Chan, *J Chrom B*, 805 (2004) 101-112.
6. M. Shipkova, E. Schutz, V. W. Armstrong, P. D. Niedmann, M. Oellerich, E. Wieland. *Clin Chem* 46:3 (2000) 365-372.



Proteos provides proteins and peptides for discovery research. We have extensive expertise in all expression systems including *E. coli*, insect cells, and mammalian cells.

What makes us unique in the industry is our commitment to:

- Your specifications
- Your timelines
- Customer communications

Deliverables include:

- Stably transfected cell lines
- 1-100mg quantities of crystallographic grade proteins
- 1-100mg of synthetic peptides

In essence, we become an extension of your bench, reproducing your protocols or conducting research to develop new protocols. We work with many clients to help successfully resolve their protein/peptide challenges by providing expert consultation with:

- Design of cDNA constructs
- Expression systems
- Purification protocols
- Custom peptides
- Assays

Let Proteos become an extension of your bench with our protein and peptide expertise!

269-372-3423 • FAX 269-353-5550 • gbgross@proteos.net • www.proteos.net

PreClinOmics, Inc

Contract Research Services and Consulting

Metabolic Disease Experts Specialized
in Providing Pre-Clinical, In-Vivo
Services to Support Drug Discovery



- Study Consulting and Management
- Protocol Development
- Model Development
- Laboratory Services

317-872-6001 www.preclinomics.com

Science education is important! The days of Bunsen burners and ring stands are fading fast.

Problem: Students are losing their curiosity about science at an early age.

Problem: Schools can't afford to give students experience with the most modern scientific instruments.

Problem: Instruments for science are developing at a rapid pace. It is hard for schools to keep up to date.

Problem: Modern instruments are not needed in schools all year long, but rather for a few weeks each year.

Problem: Maintaining instruments in good working order is not affordable for individual schools.

Solution: Science Vans

Science vans are typically college- or university-based. They allow for cost sharing among many schools by delivering well maintained modern instrumentation. They provide for centralized development and production of teaching modules that attract students and teachers alike. They encourage interactions among K-12 and college level instructors. They can help attract science majors to colleges and universities that sponsor science vans. In a knowledge-based economy, this is an opportunity not to be ignored.

Here is an example: www.chem.purdue.edu/cmobile

*Parents, teachers, government, business: **Get involved!** If your community does not have access to a science van, start asking questions, alert legislators and make noise!*



PURDUE
UNIVERSITY



This message is brought to you by the Science Education Foundation of Indiana (SEFI) and Bioanalytical Systems, Inc.